

Centri-Scope

Virtual OfficeWare enhances your GE Centricity® Physician Office - Electronic Medical Records investment with integrated value-added solutions.

Centri-Scope software accommodates physicians by allowing them to capture, identify, and embed images from a patient encounter directly into the report generated in Centricity EMR.

Real time video capture, coupled directly with the scope trigger, enables the physician to preview and obtain images during the procedure while keeping 100% hands-on the scope. Full screen, real time display allows the PC workstation to also be used as a second monitor or nurses workstation during the procedure.

Following the procedure, the physician can easily caption and resize each image. The physician is also able to select which of the captured images should be included or excluded in the final report.




Embed images into a report format and attach it to the open encounter in the EMR.

Endoscopy Procedure Image



Centri-Scope interfaces directly with the scope processor to reveal live video and to capture images.

Easily assign image positions onto the diagram.



Centri-Scope Diagram

Benefits and Value:

- Capture unlimited number of images
- Easily title images
- Indicate the position of each image on a diagram
- Select the image size from three (3) separate formats
- Attach single or multiple page reports directly to the open encounter
- Quickly produce the report including images to send to the referring physician

By using *real time* video capture, physicians are able to produce a complete account of the patient's encounter and submit a thorough report to the referring physician.